

Canadian Bioinformatics Workshops


www.bioinformatics.ca bioinformaticsdotca.github.io

Creative Commons

This page is available in the following languages: Afrikaans български Català Dansk Deutsch Ελληνικά English English (CA) English (GB) English (US) Esperanto Castellano Castellano (AR) Español (CL) Castellano (CO) Español (Ecuador) Castellano (MX) Castellano (PE) Euskara Suomeksi français français (CA) Galego אירטע hrvatski Magyar Italiano 日本語 한국어 Macedonian Metayu Nederlands Norsk Sesotho sa Leboa polski Português română slovenski jezik cpncku srpski (latiinica) Sotho svenska 中文 禅語 (台湾) isiZulu



Attribution-Share Alike 2.5 Canada




Module 2

RNA-Seq Module 2: Alignment QC

Obi Griffith and Malachi Griffith RNA-seq Analysis 2023. July 17-19, 2023

rnabio.org





Washington University in St.Louis


School of Medicine

Module 2

Learning objectives of module 3

- Visualization of RNA-seq alignments in IGV
- Alignment QC Assessment
- BAM read counting and determination of variant allele expression status

Visualization of RNA-seq alignments in IGV browser



Module 2


Alternative viewers to IGV

- Alternative viewers to IGV
 - http://www.biostars.org/p/12752/
 - http://www.biostars.org/p/71300/
- Artemis, BamView, Chipster, gbrowse2, GenoViewer, MagicViewer, Savant, Tablet, tview

Alignment QC Assessment

- 3' and 5' Bias
- Nucleotide Content
- Base/Read Quality
- PCR Artifact
- Sequencing Depth
- Base Distribution
- Insert Size Distribution

Alignment QC: 3' & 5' Bias



rnabio.org


http://rseqc.sourceforge.net/

Module 2

Gene body percentile (5'->3')

Alignment QC: Nucleotide Content

- Random primers are used to reverse transcribe RNA fragments into double-stranded complementary DNA (dscDNA)
- Causes certain patterns to be over represented at the beginning (5'end) of reads
- Deviation from expected A%=C%=G%=T%=25%



Nucleic Acids Res. 2010 Jul; 38(12): e131 Published online 2010 Apr 14. doi: 10.1093/nar/gkq224

Journal List > Nucleic Acids Res > v.38(12); 2010 Jul > PMC2896536

Biases in Illumina transcriptome sequencing caused by random hexamer

priming

Kasper D. Hansen,^{1,*} Steven E. Brenner,² and Sandrine Dudoit^{1,3}

Author information
Article notes
Copyright and License information

This article has been cited by other articles in PMC.

http://rseqc.sourceforge.net/


rnabio.org

Module 2

9


Alignment QC: Quality Distribution

- Phred quality score is widely used to characterize the quality of base-calling
- Phred quality score = -10xlog(10)P, here
 P is probability that base-calling is wrong
- Phred score of 30 means there is 1/1000 chance that the base-calling is wrong
- The quality of the bases tend to drop at the end of the read, a pattern observed in sequencing by synthesis techniques



Alignment QC: PCR Duplication


- Duplicate reads are reads that have the same start/end positions and same exact sequence
- In DNA-seq, reads/start point is used as a metric to assess PCR duplication rate
- In DNA-seq, duplicate reads are collapsed using tools such as picard
- How is RNA-seq different from DNA-seq?




Alignment QC: Sequencing Depth

• Have we sequenced deep enough?

- In DNA-seq, we can determine this by looking at the average coverage over the sequenced region. Is it above a certain threshold?
- In RNA-seq, this is a challenge due to the variability in gene abundance
- Use splice junctions detection rate as a way to identify desired sequencing depth
- Check for saturation by resampling 5%, 10%, 15%, ..., 95% of total alignments from aligned file, and then detect splice junctions from each subset and compare to reference gene model.
- This method ensures that you have sufficient coverage to perform alternative splicing analyses




Alignment QC: Base Distribution



• Your sequenced bases distribution will depend on the library preparation protocol selected

Alignment QC: Insert Size




PE reads	R1>	<r2< th=""></r2<>
fragment	~~~====================================	
insert	=======================================	
inner mate		

http://thegenomefactory.blogspot.ca/2013/08/paired-end-read-confusion-library.html

Alignment QC: Insert Size

Mean=60;SD=52



Consistent with library size selection?

BAM read counting and variant allele expression status



• A variant C->T is observed in 12 of 25 reads covering this position. Variant allele frequency (VAF) 12/25 = 48%.

• Both alleles appear to be expressed equally (not always the case) -> heterozygous, no allele specific expression

• How can we determine variant read counts, depth of coverage, and VAF without manually viewing in IGV?

We are on a Coffee Break & Networking Session

Workshop Sponsors:



Canadian Centre for Computational Genomics





g **Genome**Canada